

# **Blood RNA Miniprep Kit**

## **(BW-R6411)**

### **Contents**

Kit Contents .....	2
Introduction .....	2
Storage and Stability .....	3
Before Starting .....	3
Important Notes .....	3
Materials not Supplied .....	3
Protocol for Total RNA Extraction From White Blood Cells (Leukocytes) .....	4
Options: Removal of genomic DNA using DNase digestion .....	5
Trouble Shooting Guide .....	6
Limited Use and Warranty .....	7

## Kit Contents

Catalog#	BW-R6411-00	BW-R6411-01	BW-R6411-02
Preps	4	50	250
Buffer LY	2.4 mL	28 mL	135 mL
Buffer RB	3 mL	30 mL	135 mL
RNA Wash Buffer	2 mL	24 mL	3 x 24 mL
10 x Red Blood Cell Lysis Solution	3 mL	30 mL	135 mL
DEPC-Treated ddH <sub>2</sub> O	500 µL	10 mL	30 mL
DNA Clearance Column	4	50	250
RNA Columns	4	50	250
Collection Tubes	8	100	500
1.5 mL RNase-free microfuge tube	4	50	250
User Manual	1	1	1

\* Add 8 mL (BW-R6411-00) or 96 mL (BW-R6411-01) or 96 mL (BW-R6411-02) 96-100% ethanol to each RNA Wash Buffer bottle before use.

## Introduction

The Blood RNA Miniprep Kit provides an easy and fast method for isolating total RNA from white blood cells within 30 min. The Kit combines the reversible binding properties of ezBind RNA technology with a specially designed buffer system which can effectively remove DNA before RNA isolation. The lysate is passed through an EZgene™ DNA Clearance Column which will trap the genomic DNA. AND trace genomic DNA can be eliminated by DNase I treatment (See detail in the protocol) when it is necessary.

## Storage and Stability

All components can be stored at 4-28 °C. All kit components are guaranteed for 12 months from the date of production.

## Before Starting

Prepare all components and get all necessary materials ready by examining this user manual and become familiar with each step and pay special attention to the followings.

## Important Notes

- ☼ Add 1% volume of  $\beta$ -mercaptoethanol to Buffer LY before use and store at 4°C.
- ☼ Add 8 mL (BW-R6411-00) or 96 mL (BW-R6411-01) or 96 mL (BW-R6411-02) 100% ethanol to each RNA Wash Buffer before use.
- ☼ Red Blood Cell Lysis Solution is supplied as 10 x concentrate. Calculate the amount of buffer to be used and mix one part of 10 x Red Blood Cell Lysis Solution with 9 part of ddH<sub>2</sub>O before use.

## Materials not Supplied

- ☼ Tabletop microcentrifuge .
- ☼ Vacuum manifold if use vacuum protocol.
- ☼ 100% ethanol

***Note: Perform all steps including centrifugation at room temperature***

## Protocol for Total RNA Extraction From White Blood Cells (Leukocytes)

*Calculate the amount of 10 x Red Blood Cell Lysis Solution to be used. Mix one part of 10 x Red Blood Cell Lysis Solution with 9 part of ddH<sub>2</sub>O before use.*

1. Transfer 1 mL of whole blood (collected in heparinized or EDTA-treated tubes) into a 15 mL conical tube. Add **3 volumes** of **Red Blood Cell Lysis Solution** (see Page 3) and mix the solution by inverting the tube 5 times. Incubate on ice for 10 min.
2. Centrifuge the blood sample at 3,000 rpm (600 x g) for 5 min at 4 °C. Remove the supernatant by carefully pipetting from the top of the sample. Do not disturb the precipitated leukocyte pellets. Carefully wash the pellets with **2 mL** of **Red Blood Cell Lysis Solution** and centrifuge the sample at 3,000 rpm (600 x g) for 5 min at 4°C.
3. Remove the supernatant without disturbing the leukocyte pellets.
4. Add **500 µL Buffer LY** to the leukocytes pellets and vortex the solution for 1 min.
5. Transfer the cleared lysate to a DNA Clearance column pre-inserted in a 2 mL Collection Tube. Centrifuge at 13,000 rpm for 2 min. Discard the DNA Clearance column and save the flow-through.

***Note:** This step is for genomic DNA removal.*

6. Add **0.5 volume** 100% ethanol to the lysate (for example: 250 µL 100% ethanol for 500 µL lysate).
7. Transfer the solution into a RNA column and centrifuge at 13,000 rpm for 1 min. Discard the collection tube with the flow-through and put the column back to a new collection tube.
8. Add **500 µL Buffer RB** to the column and centrifuge at 13,000 rpm for 30s. Discard the flow-through.
9. Add **500 µL RNA Wash Buffer** to the column and centrifuge at 13,000 rpm for 1 min. Discard the flow-through.

***Ensure that ethanol has been added to RNA Wash Buffer before use.***

10. Add another **500 µL RNA Wash Buffer** to the column and centrifuge at 14,000 rpm for 30 s. Discard the flow-through.
11. Centrifuge the column at 13,000 rpm for 1 min. Discard the flow-through and collection tube.
12. Put the column into a new collection tube. Centrifuge the column, **with the lid open**, at 13,000 rpm for 2 min.

***It is critical to remove residual ethanol for optimal elution.***

13. Place the column to a RNase-free 1.5 mL tube, add **50-100 µL DEPC-treated ddH<sub>2</sub>O** to the column and centrifuge at 13,000 rpm for 2 min. The RNA is in the flow-through liquid. Store the RNA solution at -20°C.

**Note:** It is highly recommended that RNA quality be determined before downstream applications. The quality of RNA can be assessed by denatured agarose gel electrophoresis with the ethidium bromide staining. Several sharp bands should appear on the gel including 28S and 18S ribosomal RNA bands as well as certain populations of mRNA and bands. If these bands smear towards lower molecular weight RNAs, then the RNA has undergone major degradation during preparation, handling or storage, RNA molecule less than 200 bases in length do not efficiently bind to the RNA column. An  $A_{260}/A_{280}$  ratio of 1.8-2.0 corresponds to 90-100% pure nucleic acid.

## Options: Removal of genomic DNA using DNase digestion

DNA digestion is necessary for downstream applications that are sensitive to very small amounts of DNA, for example, RT-PCR with low-abundance target. Generally, it is not required to do so since the EZgene RNA purification kit selectively isolates RNA and eliminates most of the DNA. If there is DNA contamination, either reduces the tissue amount or cell .

Catalog#	BW-R6411-00	BW-R6411-01	BW-R6411-02
<b>Preps</b>	4	50	250
DNase I	25 u	260 u	1300 u
1 x DNase I Buffer	300 $\mu$ L	3 mL	15 mL
DNase Stop Buffer	200 $\mu$ L	2.4 mL	12 mL

\*DNase I , DNase Stop Buffer not supplied. They could be purchased from BEIWO

\*Add 800  $\mu$ L ( BW-R6411-00 ) or 9.6 mL (BW-R6411-01) or 48 mL (BW-R6411-02) 100% ethanol to DNase Stop Buffer before use. The final ethanol is 80% (v/v).

### Protocol for Removal of genomic DNA using DNase digestion

1. After loading the sample into the RNA column, proceed to the following step for DNase I digestion.
2. Place column in a clean 2 mL collection tube, and add **500  $\mu$ L Buffer RB**. Centrifuge as above and discard flow-through. Reuse the collection tube for next step.
3. Add **50  $\mu$ L DNase I (2U , RNase-free) Mixture** onto the middle of the column and incubate at room temperature for 15 min. Add **200  $\mu$ L DNase Stop Buffer** onto the column and centrifuge at 13,000 rpm for 1 min. Discard the flow- through. Add **300  $\mu$ L RNA Wash Buffer** to the column and centrifuge at 13,000 rpm for 1 min. Discard the flow-through.

## Trouble Shooting Guide

Problem	Possible reason	Suggested Improvement
Low $A_{260}/A_{280}$ ratios	Protein contamination	Do a Phenol: Chloroform extraction. Loss of total RNA (up to 40%) should be expected.
Low $A_{260}/A_{280}$ ratios	Guanidine Thiocyanate contamination	Add 2.5 volumes of ethanol and 0.1M NaCl (final concentration) to precipitate RNA. Incubate for 30 min at $-20^{\circ}\text{C}$ . Centrifuge at 10,000 g for 15 min at $4^{\circ}\text{C}$ . Resuspend the RNA pellet in DEPC-treated water.
Low Yield	RNA in sample degraded	Freeze samples immediately in liquid nitrogen and store at $-70^{\circ}\text{C}$ after collect it.
Low Yield	The binding capacity of the membrane in the spin column was exceeded	Use of too much tissue sample exceeding the binding capacity of spin column will cause the decreasing of total RNA yield.
Low Yield	Ethanol not added to buffer	Add ethanol to the RNA Wash Buffer and DNase Stop Buffer before purification.
Genomic DNA contamination	Too much total RNA sample was used in RT-PCR.	Reduce total RNA amount used in RT-PCR to 50-100 ng.
Genomic DNA contamination	The sample may contain too much genomic DNA.	<p>Reduce the amount of starting tissue in the preparation of the homogenate. Most tissues will not show a genomic DNA contamination problem at 30 mg or less per prep.</p> <p>Reduce cell numbers to <math>1-2 \times 10^6</math> or increase buffer volume and do multiple loadings to column.</p>

## Limited Use and Warranty

This product is intended for *in vitro* research use only. Not for use in human.

This product is warranted to perform as described in its labeling and in BEIWO's literature when used in accordance with instructions. No other warranties of any kind expressed or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by BEIWO. BEIWO's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of BEIWO, to replace the products. BEIWO shall have no liability for any direct, indirect, consequential, or incidental damage arising out of the use, the results of use, or the inability to use it product.

For technical support or learn more product information, please contact us or visit our website.



**Contact Us:** [400-115-2855](tel:400-115-2855)

[www.beiwobiomedical.com](http://www.beiwobiomedical.com)

**Customer Support:**

[market@beiwobiomedical.com](mailto:market@beiwobiomedical.com)

**Technical Support:**

[tech@beiwobiomedical.com](mailto:tech@beiwobiomedical.com)