## Kit Contents

Catalog#	DC3531-00	DC3531-01	DC3531-02
Preps	2	10	25
ezBind™ Maxi Columns	2	10	25
Buffer GC	100mL	500 mL	1250 mL
DNA Wash Buffer*	15 mL	54 mL	2×54 mL
Elution Buffer	6 mL	30 mL	60 mL
50mL Collection Tubes	2	10	25
User Manual	1	1	1

\*Add 60 mL (DC3531-00) or 216 mL (dc3531-01) or 216 mL (DC3531-02) 96-100% ethanol to each DNA Wash Buffer bottle before use.

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## Limited Use and Warranty

This product is warranted to perform as described in its labeling and in Biomiga's literature when used in accordance with instructions. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Biomiga. Biomiga's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of Biomiga, to replace the products, Biomiga shall have no liability for any direct, indirect, consequential, or incidental damage arising out of the use, the results of use, or the inability to use it product.

**\* FOR RESEARCH USE ONLY.**

## Trouble Shooting Guide

<b>Problems</b>	<b>Possible reasons</b>	<b>Suggested improvements</b>
Low DNA yield	<ol style="list-style-type: none"> <li>1. Not enough Buffer GC</li> <li>2. Agarose gel doesn't melt completely</li> <li>3. Reused electrophoresis buffer with increased pH.</li> <li>4. Fragment &lt; 200 bp</li> <li>5. Fragment &gt; 10 kb</li> </ol>	<ol style="list-style-type: none"> <li>1. Determine the volume of Buffer GC to be used correctly as instructed.</li> <li>2. Make sure to set the water bath to 55-60°C to allow gel to melt completely. Add more Buffer GC if necessary.</li> <li>3. Use fresh electrophoresis buffer.</li> <li>4. Add isopropanol as instructed.</li> <li>5. Incubate the column (after adding ddH<sub>2</sub>O or Elution Buffer) at 60°C for 15 min before elution.</li> </ol>
No DNA yield	Forgot to add ethanol to DNA Wash Buffer	Add absolute ethanol to DNA Wash Buffer as instructed before use.
DNA sample floats out of well while loading agarose gel	Ethanol was not completely removed from the column following wash step	After the wash step, centrifuge the empty column with the lid open at top speed for 1-3 min. Repeat once.
Column clogged	Agarose gel doesn't melt completely	Make sure to melt the gel at 55-60°C before loading the sample to DNA column.