

Ver: 2311

EndoFree Plasmid Midiprep Kit(Magnetic beads) (BW-MPD1420)

Kit Contents

Catalog#	BW-MPD1420-A00-01	BW-MPD1420-A00-02	BW-MPD1420-A00-03
Preps	10	25	50
Plasmid Beads	3x1.1 mL	8 mL	16 mL
Buffer A1	30 mL	65 mL	130 mL
Buffer B1	30 mL	65 mL	130 mL
Buffer N3	8 mL	20 mL	40 mL
Buffer RET	30 mL	65 mL	130 mL
Buffer KB	45 mL	110 mL	210 mL
DNA Wash Buffer*	15 mL	24 mL	54 mL
Endofree Elution Buffer	10 mL	25 mL	45 mL
RNase A (20 mg/mL)	160 µL	400 µL	760 µL
User Manual	1	1	1

*Add 60 mL(BW-MPD1420-A00-01) or 96 mL(BW-MPD1420-A00-02) or 216mL(BW-MPD1420-A00-03) 96-100% ethanol to each DNA Wash Buffer bottle before use.

Introduction

Key to the plasmid purification kit is our proprietary DNA binding system that allows the high efficient binding of DNA to our Magnetic beads while proteins and other impurities are removed by wash buffer. Nucleic acids are easily eluted with sterile water or EndoFree Elution Buffer. The purified DNA is guanidine/anion exchange resin residues free.

The EndoFree system uses a specially formulated buffer that extracts the endotoxin from the bacterial lysate. The endotoxin level is 1 to 10 EU (Endotoxin) per µg.

This kit is designed for fast and efficient purification of plasmid DNA from 15 to 50 mL of *E. coli* culture. The Magnetic beads has a DNA binding capacity of 250 µg. The purified endofree DNA is ready for downstream applications such as transfection of endotoxin-sensitive cell lines, primary cultured cells or microinjection.

The kit can be matched with a variety of automated nucleic acid extractors, such as Allsheng Auto-Pure 24A.

Important Information

Plasmid Copy Numbers: The yield of plasmid DNA depends on the origin of the replication and the size of the plasmid. The protocols are optimized for high copy number plasmid purification. For low copy number plasmids, both the culture volume and the buffer volume need to be scaled up accordingly. Please contact our customer service for further information and refer to Table 1 for the commonly used plasmids.

Table 1 Commonly used plasmids and expected yield.

Plasmid	Origin	Copy Numbers	Expected Yield (µg per 1 mL)
pSC101	pSC101	5	0.1-0.2
pACYC	p15A	10-12	0.4-0.6
pSuperCos	pMB1	10-20	0.4-1
pBR322	pMB1	15-20	0.6-1
pGEM ^R	Muted pMB1	300-400	6-7
pBluescript ^R	ColE1	300-500	6-8
pUC	Muted pMB1	500-700	8-12

Host Strains: The strains used for propagating plasmid have significant influence on yield. The strains used for propagating plasmid have significant influence on yield. Host strains such as TOP10, DH5α and C600 yield high-quality plasmid DNA. *EndA*⁺ strains such as JM101, JM110, HB101, TG1 and their derivatives, normally have low plasmid yield due to either endogenous endonucleases or high carbohydrates released during lysis. We recommended transform plasmid to an *endA*⁻ strain if the yield is not satisfactory. Please refer to Table 2 for the *endA* information.

Table 2 *endA* strains of *E. coli*.

<i>End A</i> ⁻ Strains of <i>E. coli</i>							
DH5α	DH1	DH21	JM106	JM109	SK2267	SRB	XLO
TOP10	DH10B	JM103	JM107	SK1590	MM294	Stbl2 TM	XL1-Blue
BJ5182	DH20	JM105	JM108	SK1592	Select96 TM	Stbl4 TM	XL10-Gold
<i>End A</i> ⁺ Strains of <i>E. coli</i>							
C600	JM110	RR1	ABLE [®] C	CJ236	KW251	P2392	BL21(DE3)
HB101	TG1	TB1	ABLE [®] K	DH12S TM	LE392	PR700	BL21(DE3) pLysS
JM101	JM83	TKB1	HMS174	ES1301	M1061	Q358	BMH71-18
All NM Strains				All Y Strains			

Optimal Cell Mass (OD₆₀₀ × mL of Culture): This procedure is designed for isolating plasmid

grown in standard LB medium (Luria Bertani) for 12-16 hours to a density of OD₆₀₀ 2.0 to 3.0. If rich medium such as TB or 2 × YT are used, make sure the cell density doesn't exceed 3.0 (OD₆₀₀). A high ratio of biomass over lysis buffers result in low DNA yield and purity. The Mini Column has an optimal biomass of 10-15. For example, if the OD₆₀₀ is 2.5, the optimal culture volume should be 1-5 mL. For over amount of cell numbers, either reduce the biomass or scale up the volume of Buffer A1, B1, N3 and RET.

Culture Volume: Use a flask or tube 4 times bigger in volume than the culture medium to secure optimal condition for bacteria growth. Don't exceed the maximum culture volume suggested in the protocol. Incomplete lysis due to over amount of bacterial culture results in lower yield and purity.

Storage and Stability

Buffer A1 should be stored at 4°C once RNase A is added. RNase A can be stored at room temperature (15-25 °C) and long-term storage at 4 °C. The Plasmid Beads can be transported at room temperature and stored at 4 °C for a long time. All other materials can be stored at room temperature (15-25°C). The guaranteed shelf life is 12 months from the date of production.

Before Starting

Prepare all components and get all necessary materials ready by examining this user manual and become familiar with each step and pay special attention to the followings.

Important Notes

- ☼ RNase A: 20 mg/mL. It is stable for one year at room temperature (15-25°C). Spin down RNase A vial briefly before adding to Buffer A1.
- ☼ Buffer A1 should be stored at 4°C once RNase A is added.
- ☼ Add 60 mL(BW-MPD1420-A00-01) or 96 mL(BW-MPD1420-A00-02) or 216 mL(BW-MPD1420-A00-03) 96-100% ethanol to each DNA Wash Buffer bottle before use.
- ☼ Buffer B1 precipitates below room temperature. It is critical to warm up the buffer at 37°C to dissolve the precipitates before use. Keep the cap tightly closed for Buffer B1 after use.
- ☼ Plasmid beads: Vortex well before use. It is recommended to divide according to its own use, avoid repeated opening and vortexing to reduce the magnetic bead magnetism, and increase the

magnetic bead fragments

- ☛ Ensure the availability of centrifuge capable of 12,000 ×g.
- ☛ *Carry out all centrifugations at room temperature.*

Materials not Supplied

- ☛ High speed centrifuge.
- ☛ Isopropanol
- ☛ 96-100% ethanol.
- ☛ 15 and 1.5 mL centrifugal tubes.

Safety Information

Buffer N3 contains acidic acid, wear gloves and protective eyewear when handling.

Buffer RET contains chaotropic salts, which may form reactive compounds when combines with bleach. Do not add bleach or acidic solutions directly to the preparation waste.

EndoFree Plasmid Miniprep (Magnetic beads) Protocol

1. Inoculate **15-50 mL** LB containing appropriate antibiotic with 100 μ L fresh starter culture.

Incubate at 37°C for 14-16 hours with vigorous shaking.

Note: Prolonged incubation (>16 hours) is not recommended since the *E. coli* starts to lyse and the plasmid yield may be reduced.

Note: Do not grow the culture directly from the glycerol stock.

Note: Do not use a starter culture that has been stored at 4°C.

Note: This protocol is optimized for *E. coli* strain cultured in LB medium. When using TB or 2 \times YT medium, special care needs to be taken to ensure the cell density doesn't exceed 3.0 (OD₆₀₀). Buffers need to be scaled up proportionally if over amount of cultures are being processed.

2. Harvest the bacterial culture by centrifugation for 1 minutes at 12000 \times g. Pour off the supernatant and blot the inverted tube on a paper towel to remove residue medium. Remove the residue medium completely.

Note: Residue medium will cause poor cell lysis and thus lower DNA yield.

3. Add **2.5 mL Buffer A1** (Add *RNase A* to Buffer A1 before use) and completely resuspend bacterial pellet by vortexing or pipetting.

Note: Complete resuspension is critical for bacterial lysis and lysate neutralization.

4. Add **2.5 mL Buffer B1**, mix gently by inverting the tube 10 times (*do not vortex*), and incubate at room temperature for 5 minutes.

Note: Do not incubate for more than 5 minutes.

Note: Buffer B1 precipitates (cloudy look) below room temperature. Warm up Buffer B1 at 37°C to dissolve the precipitations before use.

5. Add **800 μ L Buffer N3**, mix completely by inverting/shaking the vial for 5-10 times.

Note: Incubating the lysate in ice for 1 minute will improve the yield.

Note: It is critical to mix the solution well. If the mixture still appears conglobated, brownish or viscous, more mixing is required to completely neutralize the solution.

6. Transfer the lysate to a high speed centrifuge tube and centrifuge at 12,000 \times g for 10 minutes at room temperature. Transfer the cleared lysate to a 15 mL centrifugal tube (avoid the floating precipitates).

Note: If the rotor is cold, incubate the lysate at room temperature for 10 minutes and then perform centrifugation as described.

7. Add **2.5 mL Buffer RET** and **2.5 mL** isopropanol. Mix immediately by sharp shaking.
 8. Immediately add **300 µl plasmid beads** to the step7 **Collection Tube**. Vortex for 1 min at RT, incubate at room temperature for 5 min, then place the tube on an Magnet and allow beads to separate with the tube on the Magnet, and then discard the supernatant from the tube.
 9. Add **4 mL Buffer KB** into the **tube**, vortex for 30 s, place the tube on an Magnet and allow beads to separate with the tube on the Magnet, discard the supernatant . Put the tube back to the Magnet.
- Note:** This step is NOT necessary if the plasmid is being purified from *endA*- strain such as DH5α and TOP10. Buffer KB is necessary for *endA*+ strains such as HB101, JM110, JM101 and their derived strains.
10. Add **4 mL DNA Wash Buffer** (*Add ethanol to DNA Wash Buffer before use*) into the **tube**, vortex for 30 s, place the tube on an Magnet and allow beads to separate with the tube on the Magnet, discard the supernatant . Put the tube back to the Magnet.
 11. Add **4 mL** absolute ethanol into the **tube**, vortex for 30 s, place the tube on an Magnet and allow beads to separate with the tube on the Magnet, discard the supernatant . Put the tube back to the Magnet.
 12. Dry the beads at RT. Tube can also be dried at room temperature for 5~10 minutes.
 13. Add **500~800 µl** pre-heat(65°C) **Endofree Elution Buffer** to the tube, vortex for 1 min, and then incubate for 5 minutes at room temperature. During incubation, mix 1~2 times.
 14. Place the tube on an Magnet and allow beads to separate for 2 minutes. transfer the supernatant into a new 1.5 mL centrifugal tube.
 15. The DNA concentration can be calculated as follows,

Concentration (µg/mL)=OD₂₆₀×50×dilution factor.

Kit contents (24 well preloaded plates)

Catalog#	BW-MPD1420-A24-10	BW-MPD1420-A24-11	BW-MPD1420-A24-12	BW-MPD1420-A24-13
Preps	24 x 1	24 x 4 (96)	24 x 10 (240)	24 x 20 (480)
24 Well tip comb	1	4	10	20
Plasmid beads	300 uL	300 uL	300 uL	300 uL
Buffer A1	65 mL	245 mL	610 mL	410 mL x 3
Buffer B1	65 mL	245 mL	610 mL	410 mL x 3
Buffer N3	20 mL	70 mL	170 mL	340 mL
Buffer RET	2.5 mL	2.5 mL	2.5 mL	2.5 mL
Buffer KB	4 mL	4 mL	4 mL	4 mL
DNA Wash Buffer	4 mL	4 mL	4 mL	4 mL
RNase A (20 mg/mL)	390 uL	2 x 740 uL	4 x 920 uL	6 x 1.3 mL
Endofree Elution Buffer	500 uL	500 uL	500 uL	500 uL
User Manual	1	1	1	1

Applications: Allsheng Auto-Pure 24A

1. Take 15-50 mL of bacterial culture and follow the previous steps 1-6 to obtain the lysate supernatant.
2. Take a 24-well deep-well plate and add samples and reagents to the plate according to Table 1 below. If the kit is not pre-installed, the following reagents need to be added by yourself. (The total volume of each well must not exceed 10 mL.)

Table 1 plate set-up

24-well plate No.	Plate position	Sample / reagent	Vol. (μL)	Kit description	Note
Binding	1	lysate supernatant	5000	Added by user	Try to take the clear lysate, and then add buffer RET and isopropanol in turn to get the mixture.
		isopropanol	2500	Added by user	
		Buffer RET	2500	The reagent has been added, no need for user to add	
Beads	3	Plasmid-L Beads	300	The reagent has been added, no need for user to add	/
		ddH ₂ O/Elution buffer	700		
Wash 1	4	Buffer KB	4000	The reagent has been added, no need for user to add	/
Wash 2	5	DNA Wash Buffer	4000	The reagent has been added, no need for user to add	/



Wash 3	6	absolute ethanol	4000	Added by user	/
		24 Well tip comb	/	/	
Elution	8	Endofree Elution Buffer	500	The reagent has been added, no need for user to add	Elution volume can be adjusted according to specific requirements, at least 60ul.

- Start the instrument, place a new clean 24 well tip comb in the absolute ethanol plate (plate 6), and put 24-well plates into the corresponding position in the instrument.
- Use the program (table 2).
- Collect products after the program is completed. Take out 24-well plate, and pipette the product in 6/12 rows of Elution wells into sterile EP tube, and store at -20℃ or -80℃.

Table 2. Extraction procedures

step	name	Plate	Mix time(min)	Min Amp(%)	Wait time(min)	Vol. (uL)	Mix Speed (1-10)	Temp.(℃)	segment (0-5)	Cycle (1-10)	Mag. speed (1-10)	1 st Magnetize time (s)	2 nd Magnetize time (s)	3 rd Magnetize time (s)	Lip-1 vl (0-30 s)
1	Load	6	-	-	-	-	-	-	-	-	-	-	-	-	-
2	Mix	1	0.2	80	0	9000	7	OFF	0						
3	Beads	3	0	80	0	1000	5	OFF	1	2	1	30			
4	Binding	1	10	80	0	9000	6	OFF	3	2	1	20	30	20	0
5	Wash1	4	0.5	80	0	4000	7	OFF	3	2	1	15	20	15	0
6	Wash2	5	0.5	80	0	4000	7	OFF	3	2	1	15	20	15	0
7	Wash3	6	0.5	80	3.0	4000	7	OFF	3	2	1	15	20	15	0
8	Elute	8	5	80	0	500	2	65	2	2	1	30	30	-	30
9	Drop	3	0.2	-	-	1000	3	-	-	-	-	-	-	-	-
10	Unload	6													

Note: Set to heat up and then action, cooling action synchronization.

Purification of Low-Copy-Number Plasmid/Cosmid

The yield of low copy number plasmid is normally around 0.1-1 µg/mL of overnight culture. For isolating low copy number or medium copy number plasmid DNA, use the following guideline:

- ☛ Culture volume: Use **2×volumes** of the high copy number culture.
- ☛ Use **2×volumes** of the **Buffer A1, Buffer B1, Buffer N3 and Buffer RET**. Additional buffers can be purchased from Biomiga.
- ☛ Use **same volume** of **DNA Wash Buffer** and **EndoFree Elution Buffer**.

Trouble Shooting Guide

Problems	Possible Reasons	Suggested Improvements
Low yield	Poor cell lysis.	Resuspend pellet thoroughly by vortexing and pipetting prior to adding Buffer B1.
		Make fresh Buffer B1 if the cap had not been closed tightly. (Buffer B1 : 0.2 M NaOH and 1% SDS).
Low yield	Bacterial culture overgrown or not fresh.	Grow bacterial 12-16 hours. Spin down cultures and store the pellet at -20°C if the culture is not purified the same day. Do not store culture at 4°C overnight.
Low yield	Low copy number plasmid.	Increase culture volume and the volume of Buffer A1, B1, N3 and RET as instructed on page 9.
No DNA	Plasmid lost in host <i>E. coli</i> .	Prepare fresh culture.
Genomic DNA contamination	Over-time incubation after adding Buffer B1.	Do not vortex or mix aggressively after adding Buffer B1. Do not incubate more than 5 minutes after adding Buffer B1.
RNA contamination	RNase A not added to Buffer A1.	Add RNase A to Buffer A1.
Plasmid DNA floats out of wells while running in agarose gel	Ethanol traces were not completely removed from column.	Make sure that no ethanol residue remains in the silicon membrane before elute the plasmid DNA. Re-centrifuge or vacuum again if necessary.

Limited Use and Warranty

This product is warranted to perform as described in its labeling and in BEIWO's literature when used in accordance with instructions. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by BEIWO. BEIWO's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of BEIWO, to replace the products, BEIWO shall have no liability for any direct, indirect, consequential, or incidental damage arising out of the use, the results of use, or the inability to use it product.

For technical support or learn more product information, please contact us at 400-115-2855 or visit our website at www.beiwobiomedical.com