

EZgene™ Express Plasmid Maxiprep Kit

(BW-PD1568)

contents

Kit Contents	2
Introduction	2
Storage and Stability	3
Before Starting	3
Important	3
Important Notes	4
Materials supplied by users	5
Safety Information	6
Protocol	7
Express Plasmid Midiprep Flow Chart	10
Trouble Shooting Guide	11
Limited Use and Warranty	13

Kit Contents

Catalog#	BW-PD1568-00	BW-PD1568-01	BW-PD1568-02
Preps	2	10	25
ezBind™ Columns	2	10	25
Filter Syringe	2	10	25
2.0 ml Microfuge tube	4	20	50
Plastic wrench	1	1	1
Buffer GBL	8 mL	30 mL	70 mL
Buffer A1	22 mL	110 mL	270 mL
Buffer B1	22 mL	110 mL	270 mL
Buffer C1	27 mL	135 mL	2×170mL
DNA Wash Buffer*	15 mL	54 mL	2×54mL
RNase A (20 mg/mL)	110 µL	550 µL	1.35 mL
Elution Buffer	6 mL	30 mL	75mL
User Manual	1	1	1

*Add 60 mL (PD1568-00) or 216 mL (PD1568-01) or 216 mL (PD1568-02) 96-100% ethanol to each DNA Wash Buffer bottle before use.

Introduction

Key to the kit is our proprietary DNA binding systems that allow the highly efficient binding of DNA to our ezBind™ matrix while proteins and other impurities are removed by Wash Buffer. Nucleic acids are easily eluted with sterile water or Elution Buffer.

Unlike other kits in the markets, no chaotropic salts are contained in the buffer of our patented plasmid purification kit. The purified DNA is guanidine/anion exchange resin residues free.

This kit is designed for fast and efficient purification of plasmid DNA from 100 to 200 mL of E. coli culture. The column has a plasmid DNA binding capacity of 1000 µg.

The purified DNA is ready for high performance of downstream applications such as transfection of robust cells such as HEK293, restriction mapping, library screening, sequencing, as well as gene therapy and genetic vaccinations.

Storage and Stability

Buffer A1 should be stored at 4°C once RNase A is added. All other materials can be stored at room temperature(4-28°C). The Guaranteed shelf life is 12 months from the date of production.

Before Starting

Prepare all components and get all necessary materials ready by examining this instruction booklet and become familiar with each step.

Important

- ☼ RNase A: It is stable for more than half a year when stored at room temperature. Spin down the RNase A vial briefly. Add the RNase A solution to Buffer A1 and mix well before use. Store at 4°C.
- ☼ **Buffer B1 precipitates below room temperature. It is critical to warm up the buffer at 50°C to dissolve the precipitates before use.**
- ☼ **Buffer C1 may form precipitates below 10°C, warm up at 37°C to dissolve the precipitates before use .**
- ☼ **Keep the cap tightly closed for Buffer B1 after use.**
- ☼ Make sure the availability of centrifuge (13,000 rpm). Especially, after mixing the lysate with ethanol, the sample needs to be processed immediately by centrifugation.
- ☼ Please make sure the DNA column is tighten (using the plastic wrench) ,w hen add the lysate to DNA column.

Carry out all centrifugations at room temperature.

For certain reasons,Buffer KB are not involved in this kit from now on

Important Notes

Plasmid Copy Numbers: The yield of plasmid DNA depends on the origin of the replication and the size of the plasmid. The protocols are optimized for high copy number plasmid purification. For low copy number plasmids, both the culture volume and the buffer volume need to be scaled up 2 times. Please contact our customer service for further information and reference Table 1 for the commonly used plasmids, **Table 1 Commonly used plasmid and expected yield.**

Plasmid	Origin	Copy Numbers	Expected Yield (µg per 1 mL)
pSC101	pSC101	5	5
pACYC	P15A	10-12	5-10
pSuperCos	pMB1	10-20	10-20
pBR322	pMB1	15-20	10-20
pGEM ^R	Muted pMB1	300-400	100-150
pBluescript ^R	ColE1	300-500	100-200
pUC	Muted pMB1	500-700	150-250

Host Strains: The strains used for propagating plasmid have significant influence on yield. Host strains such as Top 10 and DH5a yield high-quality plasmid DNA. *endA*⁺ strains such as JM101, JM110, HB101, TG1 and their derivatives, normally have low plasmid yield due to either endogenous endonucleases or highcarbohydrates released during lysis. We recommend transform plasmid to an *endA*⁻ strain if the yield is not satisfactory. For purifying plasmid DNA from *endA*⁺ strains (Table 2), we recommend use product PD1712.

Table2 endA strains of E. Coli.

<i>End A- Strains of E.coli</i>							
DH5 α	DH 1	DH 21	JM1 06	JM109	SK226 7	SRB	XLO
TOP 10	DH 10B	JM 103	JM1 07	SK159 0	MM29 4	Stbl2 ™	XL1- Blue
BJ51 82	DH 20	JM 105	JM1 08	SK159 2	Select9 6™	Stbl4 ™	XL10- Gold
<i>End A+ Strains of E. coli</i>							
C600	JM1 10	RR 1	ABLE ® C	CJ236	KW251	P2392	BL21(D E3)
HB1 01	TG1	TB 1	ABLE ® K	DH12S ™	LE392	PR70 0	BL21(D E3) pLysS
JM10 1	JM8 3	TK B1	HMS1 74	ES1301	M1061	Q358	BMH 71-18
All NM strains				All Y strains			

Optimal Cell Mass (OD₆₀₀ x mL of Culture): This procedure is designed for isolating plasmid grown in standard LB medium (Luria Bertani) for 12 -16 hours to a density of OD₆₀₀ 2.0 to 3.0. If rich mediums such as TB or 2xYT are used, make sure the cell density doesn't exceed 3.0 (OD₆₀₀). A high ratio of biomass over lysis buffers result in low DNA yield and purity. The column has an optimal biomass of 150-250. For example, if the OD₆₀₀ is 3.0, the optimal culture volume should be 50 to 80 mL.

Culture Volume: Use a flask or tube 4 times bigger in volume than the culture medium to secure optimal condition for bacteria growth. Don't exceed the maximum culture volume suggested in the protocol. Incomplete lysis due to over amount of bacterial culture results in lower yield and less purity.

Materials supplied by users

☼ 100% ethanol

- ☼ High speed centrifuge
- ☼ 50 mL high speed centrifuge tubes
- ☼ 50 mL conical tubes
- ☼ Isopropanol if precipitate the plasmid DNA.

Safety Information

Buffer C1 contains acetic acid, wear gloves and protective eyewear when handling.

Buffer C1 contains chaotropic salts, which may form reactive compounds when combines with bleach. Do not add bleach or acidic solutions directly to the preparation waste.

倍乐医药

Protocol

1. **Inoculate 100-200 mL LB containing appropriate antibiotic with 100 µl fresh starter culture. Grow at 37 °C for 14-16 h with vigorous shaking.**

Note: The best way to prepare a starter culture: Inoculate a single colony from a freshly grown selective plate into 1 ml LB medium containing the appropriate antibiotic and grow at 37°C for 6 -8 h with vigorous shaking (~250 rpm). Warning: Do not use more than 100 ml culture. Need to scale up buffers if processing more than 100 mL culture.

2. **Column equilibration:** The end assembly was separated from the ezBind™ Column using a Plastic wrench and inserted into a 2 mL Microfuge Tube, and add **1 mL Buffer GBL** to ezBind™ Column. Centrifuge for 2 min at 8,000 rpm in a table-top centrifuge. Discard the tube, and reassembling the end assembly back together using the Plastic wrench. (Please use freshly treated spin column).

3. **Harvest the bacterial by centrifugation at 5,000 g for 10 min at room temperature. Pour off the supernatant and blot the inverted tube on paper towels to remove residual medium.**

Note: Complete removal of residue medium is critical for bacteria lysis in the next step.

4. **Add 10 mL Buffer A1 and completely resuspend bacterial pellet by vortexing or pipetting.** Ensure that RNase A has been added into Buffer A1 before use.

Note: Complete resuspension is critical for optimal yield.

5. **Add 10 mL Buffer B1, mix thoroughly by inverting 10 times with gentle shaking. Incubate for 5 -10 min to obtain a slightly clear lysate.** Complete lysis is critical for DNA yield. The mixture of completely lysed bacteria looks transparent.

Attention: Buffer B1 forms precipitation below room temperature, if solution becomes cloudy, warm up at 37°C to dissolve before use.

6. **Add 12 mL Buffer C1, mix completely by inverting the tube 10 times and shaking for 5 times.** It is critical to mix the solution well, if the mixture still appears conglobated, brownish or viscous; more mix is required to completely neutralize the solution.

7. **Two options for clearing the lysates:**

High speed centrifuge: Transfer the lysate to a high speed centrifuge tube and centrifuge at 13,000 rpm (14,000-18,000 g) for 15 min at room temperature! Transfer the cleared lysate to a 50 ml conical tube.

ezFilter syringe: Pour the lysate into the barrel of the filter syringe and set the syringe in a 50 mL conical tube. Incubate at room temperature for 10 min. The white precipitates should float to the top. Gently insert the plunger to expel the cleared lysate to the tube, stop when feel strong resistance, some of the lysate may remain in the flocculent precipitate.

Note: To avoid clog of the syringe: Use less than 100 mL of overnight culture and mix the lysate well after adding Buffer C1. Alternatively, transfer the lysate to another syringe filter. Additional syringe filter can be purchased from Biomiga separately.

8. Add **12 mL absolute ethanol** (96-100%) to the cleared lysate. Mix well by sharp shaking for 5 times.
9. Add the lysate/ethanol mixture into the pretreated DNA column set in a 50 mL conical tube. Use the plunger to expel the lysate through the column.

Note : Please make sure the DNA column is tightened (using the plastic wrench), when add the lysate to DNA column.

10. Disassemble the membrane column from the DNA column **using the plastic wrench**. Take the plunger out of the DNA column. Assemble the DNA column with the same membrane column and add the remaining mixture. Expel the plunger until all of mixture goes through the DNA binding membrane.

Note: Please make sure to disassemble the membrane column from the DNA column before taking the plunger out of the column. Otherwise, the pressure difference between the inside and outside of the column would break the membrane.

11. Add **10 mL DNA Wash Buffer** to the assembled DNA column (disassemble the membrane column from the DNA column before taking the plunger out), expel the buffer out with the plunger.
12. Use the plastic wrench to detach the end component from the midiprep

column and insert it into a 2.0 ml eppendorff tube.

13. Spin the column at 13,000- 15,000 rpm (Max speed) for 1 min. Decant the flow through and put the column back to the tube. Spin the column at max speed for 2 min, and then air-dry the membrane column at room temperature for 3-10min.
14. Repeat step “14”
15. Insert the membrane column into a new sterile 1.5mL or 2.0mL Microfuge tube. Add 500 μ L Elution Buffer (Prewarm the Elution Buffer 60°C increases the yield) to the center of the column and incubate for 1 min at room temperature. Elute the DNA by centrifugation at 12,000 rpm for 1 min. Then add 300 μ L Elution Buffer to the center of the column and incubate for 1 min at room temperature. Elute the DNA by centrifugation at 12,000 rpm for 1 min.
16. Optional: Add the eluted DNA back to the column for another elution. The first elution normally yields about 60-70% of the DNA and the second elution yields another 20-30% of the DNA.
Note: The DNA is ready for downstream applications such as cloning or transfection of HEK293 cells.

Note: It's highly recommended to remove the endotoxin if the DNA is used for endotoxin-sensitive cell lines, primary cultured cells or microinjection. Please use Biomiga's Endofree Plasmid Midiprep Kit (PD1420 and PD1422) for transfecting endotoxin-sensitive cell lines, primary cultured cells or microinjection.

Express Plasmid Midiprep Flow Chart



Lysis and neutralization

Clear Lysate through a
Syringe filter

Bind DNA to the DNA column

Wash the column with DNA Wash Buffer
by a plunger

Detach the column and insert to a 2.0 mL tube.
Spin the column to dry and elute the DNA.

Trouble Shooting Guide

Problems	Possible Reasons	Suggested Improvements
Low Yield	Poor Cell lysis.	<ul style="list-style-type: none"> Resuspend pellet thoroughly by vortexing and pipetting prior adding buffer B1. Make fresh buffer B1 if the cap had not been closed tightly. (Buffer B1: 0.2M NaOH and 1%SDS).
Low Yield	Bacterial culture overgrown or not fresh.	Grow bacterial 12-16 hours. Spin down cultures and store the pellet at -20°C . if the culture is not purified the same day. Do not store culture at 4°C over night.
Low Yield	Low copy-number plasmid.	Increase culture volume Increase the volume of buffer A1, B1, C1 and 100% ethanol proportionally with the ratio of 1:1:1.2:1.2.
No DNA	Plasmid lost in Host <i>E.coli</i>	Prepare fresh culture.
Genomic DNA contamination	Over-time incubation after adding buffer B1.	Do not vortex or mix aggressively after adding buffer B1. Do not incubate more than 5 minutes after adding solution B1.
RNA contamination	RNase A not added to solution A1.	Add RNase A to buffer A1.

BW-PD1568 EZgene™ Express Plasmid Maxiprep Kit

Plasmid DNA floats out of wells while running in agarose gel, DNA doesn't freeze or smell of ethanol	Ethanol traces not completely removed from column.	Make sure that no ethanol residual remaining in the silicon membrane before eluting the plasmid DNA. Re-centrifuge or vacuum again if necessary.
Low Yield	Poor Cell lysis.	<ul style="list-style-type: none"> • Resuspend pellet thoroughly by vortexing and pipetting prior adding buffer B1. • Make fresh buffer B1 if the cap had not been closed tightly. (Buffer B1: 0.2M NaOH and 1%SDS).

Limited Use and Warranty

This product is intended for *in vitro* research use only. Not for use in human.

This product is warranted to perform as described in its labeling and in BEIWO's literature when used in accordance with instructions. No other warranties of any kind expressed or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by BEIWO. BEIWO's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of BEIWO, to replace the products. BEIWO shall have no liability for any direct, indirect, consequential, or incidental damage arising out of the use, the results of use, or the inability to use it product.

For technical support or learn more product information, please contact us or visit our website.



Contact Us: 400-115-2855

www.beiwobiomedical.com

Customer Support:

market@beiwobiomedical.com

Technical Support:

tech@beiwobiomedical.com