



8201男款钢皮 / 渐进多焦点防蓝光变色镜

可戴着走路开车+变色



远近三用 智能老花镜

 (136mm)

 (54mm)

 (42mm)

 (18mm)

 (140mm)

 (29.5g)

每款老花镜均可选: +100°/+150°/+200°/+250°/+300°/+350°/+400°



8601 钢皮男款/防蓝光双光隐形10层镀膜

隐形双光

防蓝光防辐射



可带着开车走路



每款老花镜均可选：+100°/+150°/+200°/+250°/+300°/+350°/+400°



8602 半框女款/防蓝光双光隐形10层镀膜

隐形双光

防蓝光防辐射



可带着开车走路



每款老花镜均可选: +100°/+150°/+200°/+250°/+300°/+350°/+400°



8305 半框男款/防蓝光10层镀膜

防蓝光抗疲劳

防蓝光抗辐射 国标质检合格



缓解疲劳

时尚显年轻老花镜



每款老花镜均可选: +100°/+150°/+200°/+250°/+300°/+350°/+400°



8302 半框女款/防蓝光10层镀膜

防蓝光抗疲劳

防蓝光抗辐射 国标质检合格



缓解疲劳

时尚显年轻老花镜



每款老花镜均可选: +100°/+150°/+200°/+250°/+300°/+350°/+400°



8304 半框女款/防蓝光10层镀膜

防蓝光抗疲劳

防蓝光抗辐射 国标质检合格



缓解疲劳

时尚显年轻老花镜



每款老花镜均可选: +100°/+150°/+200°/+250°/+300°/+350°/+400°



8501TR 男款/负离子防蓝光10层镀膜

平均释放负氧离子
1785个/CM³



持续释放负氧离子

防蓝光老花镜



每款老花镜均可选: +100°/+150°/+200°/+250°/+300°/+350°/+400°



8502TR 女款/负离子防蓝光10层镀膜

平均释放负氧离子
1785个/CM³



持续释放负氧离子

防蓝光老花镜



每款老花镜均可选: +100°/+150°/+200°/+250°/+300°/+350°/+400°



引领高端时尚老花品牌 —— 康强

高新技术企业证书

企业名称：温州亮彩眼镜科技有限公司

证书编号：GR202133001195

发证日期：2021年12月16日

有效期：三年

批准机关：



检验报告

报告编号：BZ22030017

样品名称：老视镜

委托单位：温州亮彩眼镜科技有限公司

检验类别：委托检验

注意事项

1. 报告为“检验测试专用章”无效。
Report is invalid without the seal only for test.
2. 未经本机构批准，不得复制（全文复制或除外）报告。
No copy (except full copy) of the report is allowed without the approval of the agency.
3. 报告无编制、审核、批准人签字无效。
Report is invalid without signature.
4. 报告涂改无效。
Report is invalid if altered.
5. 本报告的检测结果仅对样品有效。
It's effect that the result of this report relate only to the sample tested.
6. 若有异议，请于收到报告之日起十五日内向本中心提出，逾期不予受理。
If you have any dissent from the report, please raise it in fifteen days to Zhejiang Inspection & Testing Centre For ophthalmic Optics products after receiving the report or we won't deal with the case.
7. 中文与英文解释有差异的，以中文为准。
If there exist any difference between Chinese and English interpretation, take the Chinese as standard.
8. 本报告仅对委托负责，不得用于广告宣传。
The report takes responsibility only for samples tested. It can not be used for advertising.

地址：浙江省温州市瓯海区经济开发区娄桥工业园奥光路2号 邮政编码：325041
Add: Jigong Road 2, Louqiao Industrial Park, Ouhai Economic Developing Zone, Wenzhou Zhejiang Post Code: 325041

电话：0577-56575903 传真：0577-56575900
Tel: Fax:
E-mail: E-mail:
Network Address: Network Address:

浙江省眼镜产品质量检验中心

Zhejiang Inspection & Testing Centre For ophthalmic Optics products

检验报告

Inspection report

报告编号：BZ22030017 Report No. 第1页 共3页 Page 1 of 3 pages

样品名称 Name of sample(s)	型号规格 Type of sample(s)	数量 Quantity	委托单位 Client	生产单位 Manufacturer	样品数量 Quantity of sample(s)	样品编号 Sample(s) No.	样品状况 Condition of sample(s)	检测日期 Date of receipt	检测日期 Date of inspection
老视镜	+1.00D, OGD-62mm	4副	温州亮彩眼镜科技有限公司	温州亮彩眼镜科技有限公司	4副	BZ22030017- (1#-4#)	完好	2022-03-25	2022-03-25-2022-03-29
检验项目 Items of inspection	视觉测试区								
检验结论 Inspection result	依据委托要求，依据GB/T 13511.3-2019《配装眼镜 第3部分：单光老视镜》标准对样品进行检测，检验项目检测结果符合标准要求。								
备注 Annotation									

批准：黄涛 审核：吴品 编制：徐科
Approved by: 黄涛 Verified by: 吴品 Typed by: 徐科

序号	检测项目	标准要求	检测结果	单项结论
1	标志	GB/T 13511.3-2019 6	符合	合格
2	镜片材料和表面质量	GB/T 13511.3-2019 4.1.1	符合	合格
3	偏光性能	GB/T 13511.3-2019 4.1.3	符合	合格
4	镜片与镜架的几何形状	GB/T 13511.3-2019 4.1.5	符合	合格
5	厚度	GB/T 13511.3-2019 4.1.1	±0.12	合格
6	内径	GB/T 13511.3-2019 4.1.1	±0.09	合格
7	光学中心水平偏差	GB/T 13511.3-2019 4.2.2	±2.0	合格
8	光学中心垂直偏差	GB/T 13511.3-2019 4.2.3	±1.0	合格
9	光学中心垂直互差	GB/T 13511.3-2019 4.2.4	≤1.0	合格

序号	检测项目	标准要求	检测结果	单项结论
10	包装结合力	GB/T 13511.3-2019 4.1.3	符合	合格
11	抗汗性能	GB/T 13511.3-2019 4.1.3	符合	合格
12	阻燃性	GB/T 13511.3-2019 4.1.1	符合	合格
13	镜片夹持力	GB/T 13511.3-2019 4.1.4	符合	合格
14	透射比要求 (%)	GB/T 13511.3-2019 4.1.2	符合	合格

浙江省眼镜产品质量检验中心
Zhejiang Inspection & Testing Centre For ophthalmic Optics products
(温州市瓯海区质量技术监督检测所)